


softFlow Clinic Applicators

Patient interface for clinical application in adult and adolescent patients














- ✓ High wearing comfort due to individual fit, simple handling and pleasant skin tolerance
- ✓ Soft and smooth thanks to medical silicone
- ✓ Patented, integrated heating up to the nose prevents condensation
- ✓ 7 different applicator types for optimal comfort at any flow rate
- ✓ Optional interface for tracheotomy patients with flow rates up to 50 l/min



soft Flow Clinic Applicators

TECHNICAL DATA:

Type:	Applicator Clinic Small	Applicator Clinic Small-Plus	Applicator Clinic Standard	Applicator Clinic Standard-Plus	Applicator Clinic Medium	Applicator Clinic Large	Applicator Clinic Tracheal Interface
Article number:	40630001	40630004	40630002	40630005	40630010	40630013	40630019
Hygiene:	Single patient multiple use	Single patient multiple use	Single patient multiple use	Single patient multiple use	Single patient multiple use	Single patient multiple use	Single patient multiple use
Maximum flow rate:	20 l/min	45 l/min	25 l/min	35 l/min	60 l/min	60 l/min	50 l/min
Operating life:	360 therapy hours	360 therapy hours	360 therapy hours	360 therapy hours	360 therapy hours	360 therapy hours	360 therapy hours
Distance NC (A): 	12,5 mm	12,5 mm	15 mm	15 mm	15 mm	16 mm	-
Length NC (B): 	11 mm	11 mm	13 mm	13 mm	13 mm	13 mm	-
Diameter NC (D): 	3,2 mm	3,2 mm	4,0 mm	4,0 mm	4,0 mm	5,0 mm	 22 mm
System connection:							
Total length (L):	225 cm	225 cm	225 cm	225 cm	225 cm	225 cm	225 cm
Supply tube length: Inner diameter: Outside diameter:	180 cm 6 mm 9 mm	180 cm 8 mm 12 mm	180 cm 6 mm 9 mm	180 cm 8 mm 12 mm	180 cm 8 mm 12 mm	180 cm 8 mm 12 mm	225 cm 7 mm 10 mm
Loop (S) length: Inner diameter: Outside diameter:	2x45 mm 3,6 mm 5 mm	2x45 mm 5 mm 7 mm	2x45 mm 3,6 mm 5 mm	2x45 mm 3,6 mm 5 mm	2x45 mm 5 mm 7 mm	2x45 mm 5 mm 7 mm	-
Heating:	max. 24 VA	max. 24 VA	max. 24 VA	max. 24 VA	max. 24 VA	max. 24 VA	max. 24 VA
CE Certificate:	CE 0297	CE 0297	CE 0297	CE 0297	CE 0297	CE 0297	CE 0297
Disposal:	Domestic waste	Domestic waste	Domestic waste	Domestic waste	Domestic waste	Domestic waste	Domestic waste
Manufacturer:	TNI medical AG Tiefe Gasse 11-13 97084 Wuerzburg Germany	TNI medical AG Tiefe Gasse 11-13 97084 Wuerzburg Germany	TNI medical AG Tiefe Gasse 11-13 97084 Wuerzburg Germany	TNI medical AG Tiefe Gasse 11-13 97084 Wuerzburg Germany	TNI medical AG Tiefe Gasse 11-13 97084 Wuerzburg Germany	TNI medical AG Tiefe Gasse 11-13 97084 Wuerzburg Germany	TNI medical AG Tiefe Gasse 11-13 97084 Wuerzburg Germany

