


# softFlow Homecare Applicators

Patient interface for homecare application in adult and adolescent patients














- ✓ High wearing comfort due to individual fit, simple handling and pleasant skin tolerance
- ✓ Soft and smooth thanks to medical silicone
- ✓ Patented, integrated heating up to the nose prevents condensation
- ✓ 7 different applicator types for optimal comfort at any flow rate
- ✓ Optional interface for tracheotomy patients with flow rates up to 50 l/min



# soft Flow Homecare Applicators

## TECHNICAL DATA:

Type:	Applicator Homecare Small	Applicator Homecare Small-Plus	Applicator Homecare Standard	Applicator Homecare Standard-Plus	Applicator Homecare Medium	Applicator Homecare Large	Applicator Homecare Tracheal Interface
<b>Article number:</b>	40630101	40630104	40630102	40630105	40630110	40630113	40630119
<b>Hygiene:</b>	Single patient multiple use	Single patient multiple use	Single patient multiple use	Single patient multiple use	Single patient multiple use	Single patient multiple use	Single patient multiple use
<b>Maximum flow rate:</b>	20 l/min	45 l/min	25 l/min	35 l/min	60 l/min	60 l/min	50 l/min
<b>Operating life:</b>	720 therapy hours	720 therapy hours	720 therapy hours	720 therapy hours	720 therapy hours	720 therapy hours	720 therapy hours
<b>Distance NC (A):</b> 	12,5 mm	12,5 mm	15 mm	15 mm	15 mm	16 mm	-
<b>Length NC (B):</b> 	11 mm	11 mm	13 mm	13 mm	13 mm	13 mm	-
<b>Diameter NC (D):</b> 	3,2 mm	3,2 mm	4,0 mm	4,0 mm	4,0 mm	5,0 mm	 22 mm
<b>System connection:</b>							
<b>Total length (L):</b>	225 cm	225 cm	225 cm	225 cm	225 cm	225 cm	225 cm
<b>Supply tube length:</b> <b>Inner diameter:</b> <b>Outside diameter:</b>	180 cm 6 mm 9 mm	180 cm 8 mm 12 mm	180 cm 6 mm 9 mm	180 cm 8 mm 12 mm	180 cm 8 mm 12 mm	180 cm 8 mm 12 mm	225 cm 7 mm 10 mm
<b>Loop (S) length:</b> <b>Inner diameter:</b> <b>Outside diameter:</b>	2x45 mm 3,6 mm 5 mm	2x45 mm 5 mm 7 mm	2x45 mm 3,6 mm 5 mm	2x45 mm 3,6 mm 5 mm	2x45 mm 5 mm 7 mm	2x45 mm 5 mm 7 mm	-
<b>Heating:</b>	max. 24 VA	max. 24 VA	max. 24 VA	max. 24 VA	max. 24 VA	max. 24 VA	max. 24 VA
<b>CE Certificate:</b>	CE 0297	CE 0297	CE 0297	CE 0297	CE 0297	CE 0297	CE 0297
<b>Disposal:</b>	Domestic waste	Domestic waste	Domestic waste	Domestic waste	Domestic waste	Domestic waste	Domestic waste
<b>Expected service life:</b>	see label	see label	see label	see label	see label	see label	see label
<b>Manufacturer:</b>	TNI medical AG Tiefe Gasse 11-13 97084 Wuerzburg Germany	TNI medical AG Tiefe Gasse 11-13 97084 Wuerzburg Germany	TNI medical AG Tiefe Gasse 11-13 97084 Wuerzburg Germany	TNI medical AG Tiefe Gasse 11-13 97084 Wuerzburg Germany	TNI medical AG Tiefe Gasse 11-13 97084 Wuerzburg Germany	TNI medical AG Tiefe Gasse 11-13 97084 Wuerzburg Germany	TNI medical AG Tiefe Gasse 11-13 97084 Wuerzburg Germany

